



## CHANGES TO HOPAID® AND ITS PACKAGING

### Why have we changed our packaging?

We will be changing the way we package Hopaid for two reasons. First, feedback from our customers and second, a necessary change in pasteurisation suppliers.

### What packaging will we be using in the future?

We are moving to a bag-in-box configuration as shown below. The new packaging includes a tap attached to a collapsible bag ensuring there is no air ingress as Hopaid is dosed. Sizes available will be 2.25 and 10 Kg. Shelf life will remain as 2 years and open shelf life will be set at 4 weeks (instead of "Use immediately once open") and improvement possible due to the new configuration.



Figure 1. Hopaid bag in box

### What are the advantages of the new presentation?

- Easier and more hygienic dosing.
- Tap can be connected to an adaptor allowing for easier process automation. See a video here: <https://www.youtube.com/watch?v=W6vKGS1v5tw>

**Note: tap adaptor is not provided as a standard.**



### Where do I get more information?

You can contact our Sales or Quality Department for more information on the new packaging.